



HyPort Series

Medical supply unit

Care for health with unrivaled flexibility



Sydney | Melbourne | Brisbane | Perth


☎ 1300HPAUST
✉ info@hpaust.com
🌐 www.hpaust.com



MRB066V123

Medical Supply Unit

With the continuous improvement of medical levels, people pay more attention to efficient clinical workflow. Mindray innovation integrates high-end technology with ergonomic design.It develops a new generation of medical supply unit HyPort series, providing an efficient, easy-to-use and safe operating room environment.



HyPort P30
Mechanical Pendant



HyPort P60
Motorized Pendant



HyPort P90
Vertically Lifting Pendant

Customized Design with Future-proof Flexibility



Comfortable layout



Easy to use



Improve efficiency

- Various supporting solutions towards different medical devices
- Innovative cable management provides a clean and orderly workstation
- Ergonomic functionality with optimized workflow efficiency
- Infection control and safety system secure the surgical progress

Safety

Safety is the first priority for medical devices, especially pendants that lift and support other heavy machines.



Brake system
Besides the standard friction brake, pneumatic or electromagnetic can be an option for additional safety.



Oxygen concentration monitoring
24-hour continuous monitoring of oxygen concentration with built-in alarm system



Infection control

- Screwless design ensures no snags or cuts
- Antimicrobial powder coating surface treatment

Anesthesia Pendant Solution

Regardless of the adoption of anesthesia solutions, it gives physicians a reliable platform to work with anesthesia machines and creates a comfortable, safe operating room environment.



With a motorized arm/lifting system, HyPort P90 can lock your anesthesia workstation or endoscopy cart and lift them off the floor. The devices can be conveniently moved into the desired area and securely parked, providing a clean and efficient workflow.



Mindray also provides a compact anesthesia pendant design, namely P30/P60 with head distribution only, ensuring a large operating space for the surgical team.



Endoscope Pendant Solution

Minimally invasive surgery is known as a technique with fewer complications, shorter hospital stays and less pain. Cameras, computers, and instruments are capable of providing very precise movements to endure the procedures. Mindray endoscope pendant with the new Σ distribution is an innovative solution for those numerous devices and cables.

Endoscope Pendant Solution

Cable management



Electrical cables hidden inside the column cavity



Built-in cable management with flexible position

Smart storage



A better basket to hold foot switches with cables




Dedicated holder for laparoscopic camera

Future-proof flexibility




Isolated cavity design for safety concern



Extendable shelf design provides more space for large equipment.


Surgical Pendant Solution

The design of the column will maximize the configuration requirements towards the quantity of electrical sockets and gas outlets. In addition, the front and rear panels can be expanded as the component panels.



Imaging Mounting Solution

Due to different clinical requirements and various display selections, the manufacturers were required to do lots of time-consuming and expensive customization previously. The Mindray imaging mounting solution provides flexible configurations.



The monitor hanger can be adjusted according to the size of the monitor to meet different viewing needs, compatible with large format screens up to 65".



Integrated information solution for Hybrid OR, Catheterization Lab, MIS OR, four screens up to 32", six screens up to 27".